

# ZNF639 antibody

Catalog Number: 117014

## Product name

ZNF639 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF639 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of ZNF639 recombinant protein (Accession Number: XM\_017006551). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117014(ZNF639 Antibody) at dilution of 1:300



IP Result of anti-ZNF639 (IP:Catalog No:117014, 4ug; Detection:Catalog No:117014 1:300) with HeLa cells lysate 3960ug.