

ZNF645 antibody



Catalog Number: 117015

Product name

ZNF645 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF645 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of ZNF645 recombinant protein (Accession Number: NM_152577). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117015(ZNF645 Antibody) at dilution of 1:1000