

# ZNF689 antibody

Catalog Number: 117021

---

## Product name

ZNF689 antibody

## Immunogen

[Human ZNF689 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF689 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse testis tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of ZNF689 recombinant protein (Accession Number: NM\_138447). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117021(ZNF689 Antibody) at dilution of 1:600