

# ZNF703 antibody



Catalog Number: 117024

## Product name

ZNF703 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZNF703 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human brain tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 58kd

## Preparation

This antibody was obtained by immunization of ZNF703 recombinant protein (Accession Number: NM\_025069). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117024(ZNF703 antibody) at dilution of 1:500