

# ZNF71 antibody



Catalog Number: 117025

---

## Product name

ZNF71 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF71 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of ZNF71 recombinant protein (Accession Number: XM\_017027080). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117025(ZNF71 Antibody) at dilution of 1:1000