

ZNF740 antibody

Catalog Number: 117028

Product name

ZNF740 antibody

Immunogen

[Human ZNF740 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF740 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, mouse liver tissue, mouse testis tissue, PC-3 cells. Observed molecular weight by Western-blot: 22-30 kDa

Preparation

This antibody was obtained by immunization of ZNF740 recombinant protein (Accession Number: NM_001004304). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117028(ZNF740 Antibody) at dilution of 1:600



IP Result of anti-ZNF740 (IP:Catalog No:117028, 3ug; Detection:Catalog No:117028 1:600) with mouse brain tissue lysate 4000ug.