

ZNF746 antibody



Catalog Number: 117029

Product name

ZNF746 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF746 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 80-90 kDa

Preparation

This antibody was obtained by immunization of ZNF746 recombinant protein (Accession Number: NM_152557). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117029(ZNF746 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:117029(ZNF746 Antibody) at dilution of 1:50 (under 40x lens)