

ZNF774 antibody

Catalog Number: 117033

Product name

ZNF774 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF774 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 55-63 kDa

Preparation

This antibody was obtained by immunization of ZNF774 recombinant protein (Accession Number: NM_001004309). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117033(ZNF774 Antibody) at dilution of 1:600



IP Result of anti-ZNF774 (IP:Catalog No:117033, 3ug; Detection:Catalog No:117033 1:600) with HepG2 cells lysate 2400ug.