

# ZFP2 antibody



Catalog Number: 117039

## Product name

ZFP2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZFP2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Neuro-2a cells. Observed molecular weight by Western-blot: 53-55 kDa

## Preparation

This antibody was obtained by immunization of ZFP2 recombinant protein (Accession Number: NM\_030613). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with (ZFP2 Antibody) at dilution of 1:1000