

ZGPAT antibody



Catalog Number: 117051

Product name

ZGPAT antibody

Specificity

human, mouse; other species not tested.

Antibody description

ZGPAT Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HepG2 cells, mouse liver tissue.

Preparation

This antibody was obtained by immunization of ZGPAT recombinant protein (Accession Number: NM_181485). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:117051(ZGPAT antibody) at dilution of 1:1000