

ZHX2 antibody



Catalog Number: 117053

Product name

ZHX2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZHX2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human heart tissue, HeLa cells, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse ovary tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 92-100 kDa

Preparation

This antibody was obtained by immunization of ZHX2 recombinant protein (Accession Number: XM_011516931). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:117053(ZHX2 antibody) at dilution of 1:800



IP Result of anti-ZHX2 (IP:Catalog No:117053, 4ug; Detection:Catalog No:117053 1:800) with mouse brain tissue lysate 2800ug.