

ZIC3 antibody



Catalog Number: 117054

Product name

ZIC3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZIC3 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, mouse cerebellum tissue, SH-SY5Y cells, Y79 cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of ZIC3 recombinant protein (Accession Number: NM_003413). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:117054(ZIC3 Antibody) at dilution of 1:600