

# ZIC3 antibody



Catalog Number: 117054

---

## Product name

ZIC3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZIC3 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, mouse cerebellum tissue, SH-SY5Y cells, Y79 cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of ZIC3 recombinant protein (Accession Number: NM\_003413). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:117054(ZIC3 Antibody) at dilution of 1:600