

# ZIK1 antibody

Catalog Number: 117056

## Product name

ZIK1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZIK1 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50-55 kDa

## Preparation

This antibody was obtained by immunization of ZIK1 recombinant protein (Accession Number: NM\_001010879). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117056(ZIK1 Antibody) at dilution of 1:300



IP Result of anti-ZIK1 (IP:Catalog No:117056, 4ug; Detection:Catalog No:117056 1:300) with Jurkat cells lysate 2800ug.