

ZIM3 antibody

Catalog Number: 117057

Product name

ZIM3 antibody

Specificity

Human; other species not tested.

Antibody description

ZIM3 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, A549 cells. Observed molecular weight by Western-blot: 55-65 kDa

Preparation

This antibody was obtained by immunization of ZIM3 recombinant protein (Accession Number: NM_052882). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117057(ZIM3 Antibody) at dilution of 1:300