

ZKSCAN1 antibody



Catalog Number: 117064

Product name

ZKSCAN1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZKSCAN1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, mouse testis tissue, PC-3 cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of ZKSCAN1 recombinant protein (Accession Number: NM_001346581). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117064(ZKSCAN1 Antibody) at dilution of



1:600

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117064(ZKSCAN1 Antibody) at dilution of 1:600