

ZNF30 antibody

Catalog Number: 117158

Product name

ZNF30 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF30 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse spleen tissue, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 65-71 kDa

Preparation

This antibody was obtained by immunization of ZNF30 recombinant protein (Accession Number: XM_011527443). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:50-1:500

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:117158(ZNF30 Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:117158(ZNF30 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)