

# ZNF323 antibody



Catalog Number: 117161

## Product name

ZNF323 antibody

## Immunogen

[Human ZNF323 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF323 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse lung tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of ZNF323 recombinant protein (Accession Number: NM\_001243244). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117161(ZNF323 Antibody) at dilution of 1:600