

ZNF391 antibody

Catalog Number: 117170

Product name

ZNF391 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF391 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of ZNF391 recombinant protein (Accession Number: XM_006715086). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-ZNF391 (IP:Catalog No:117170, 3ug; Detection:Catalog No:117170 1:300) with mouse testis tissue lysate 4000ug.