

ZNF394 antibody



Catalog Number: 117171

Product name

ZNF394 antibody

Immunogen

[Human ZNF394 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF394 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, fetal human brain tissue, HepG2 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 50-60 kDa

Preparation

This antibody was obtained by immunization of ZNF394 recombinant protein (Accession Number: NM_032164). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117171(ZNF394 Antibody) at dilution of 1:600