

ZNF397 antibody



Catalog Number: 117173

Product name

ZNF397 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

ZNF397 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse testis tissue. Observed molecular weight by Western-blot: 60-63 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of ZNF397 recombinant protein (Accession Number: NM_032347). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117173(ZNF397 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality