

ZNF410 antibody



Catalog Number: 117174

Product name

ZNF410 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNF410 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 67 kDa and 49 kDa

Preparation

This antibody was obtained by immunization of ZNF410 recombinant protein (Accession Number: NM_001242926). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117174(ZNF410 antibody) at dilution of 1:500