

# ZNF443 antibody

Catalog Number: 117176

---

## Product name

ZNF443 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse; other species not tested.

## Applications

ELISA, WB

## Antibody description

ZNF443 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, mouse lung tissue. Observed molecular weight by Western-blot: 72kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of ZNF443 recombinant protein (Accession Number: NM\_005815). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:117176(ZNF443 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality