

ZNF460 antibody



Catalog Number: 117179

Product name

ZNF460 antibody

Immunogen

[Human ZNF460 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF460 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of ZNF460 recombinant protein (Accession Number: NM_006635). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117179(ZNF460 Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:117179(ZNF460 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)