

ZNF486 antibody



Catalog Number: 117183

Product name

ZNF486 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF486 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of ZNF486 recombinant protein (Accession Number: BC008936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117183(ZNF486 Antibody) at dilution of 1:600