

# ZNF488 antibody



Catalog Number: 117184

## Product name

ZNF488 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZNF488 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of ZNF488 recombinant protein (Accession Number: BC026895). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:117184(ZNF488 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human brain using Catalog No:117184(ZNF488 antibody) at dilution of 1:50 (under 40x lens)