

# ZNF828 antibody

Catalog Number: 117250

## Product name

ZNF828 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF828 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells.

## Preparation

This antibody was obtained by immunization of ZNF828 recombinant protein (Accession Number: NM\_032436). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:500-1:5000

## Validations



IP Result of anti-ZNF828 (IP:Catalog No:117250, 4ug; Detection:Catalog No:117250 1:1000) with Jurkat cells lysate 4000ug.