

ZNF846 antibody



Catalog Number: 117252

Product name

ZNF846 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

ZNF846 Rabbit Polyclonal antibody. Positive WB detected in C6 cells. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of ZNF846 recombinant protein (Accession Number: BC118657). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



C6 cells were subjected to SDS PAGE followed by western blot with Catalog No:117252(ZNF846 Antibody) at dilution of 1:600