

# ZnT3 antibody

Catalog Number: 117256

## Product name

ZnT3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZnT3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 42kd, 48kd

## Preparation

This antibody was obtained by immunization of ZnT3 recombinant protein (Accession Number: NM\_003459). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117256(SLC30A3 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:117256(SLC30A3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:117256(SLC30A3 antibody) at dilution of 1:100 (under 40x lens)