

Znt5 antibody

Catalog Number: 117257

Product name

Znt5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Znt5 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 84-87 kDa

Preparation

This antibody was obtained by immunization of Znt5 recombinant protein (Accession Number: NM_024055). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-SLC30A5 (IP:Catalog No:117257, 4ug; Detection:Catalog No:117257 1:500) with mouse liver tissue lysate 4000ug.