

ZO-1 antibody



Catalog Number: 117260

Product name

ZO-1 antibody

Specificity

Human; other species not tested.

Antibody description

ZO-1 Rabbit Polyclonal antibody. Positive IF detected in HUVEC cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003257).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of HUVEC cells using Catalog No:117260(TJP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)