

ZO-2 antibody

Catalog Number: 117261

Product name

ZO-2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZO-2 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, MCF7 cells, mouse heart tissue. Observed molecular weight by Western-blot: 150kd

Preparation

This antibody was obtained by immunization of ZO-2 recombinant protein (Accession Number: NM_004817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117261(TJP2 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells using Catalog No:117261(TJP2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)