

# ZSCAN20 antibody



Catalog Number: 117266

## Product name

ZSCAN20 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZSCAN20 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Jurkat cells. Positive IHC detected in human small intestine tissue.

Observed molecular weight by Western-blot:  
100-110 kDa

## Preparation

This antibody was obtained by immunization of ZSCAN20 recombinant protein (Accession Number: BC008827). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117266(ZSCAN20 antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:117266(ZSCAN20 Antibody) at dilution of 1:100