

Anti-APLP1 antibody



Catalog Number: 141742

Product name

Anti-APLP1 antibody

Immunogen

[Mouse APLP1 / Amyloid-like protein 1 \(His Tag\)](#)

Specificity

Mouse APLP1 / Amyloid-like protein 1

Antibody description

Rabbit Polyclonal to Mouse APLP1

Preparation

Produced in rabbits immunized with purified, recombinant Mouse APLP1 / Amyloid-like protein 1. APLP1 / Amyloid-like protein 1 specific IgG was purified by Mouse APLP1 / Amyloid-like protein 1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse APLP1 / Amyloid-like protein 1.

IHC-P: 0.01-1µg/mL

Validations



Mouse APLP1 / Amyloid-like protein 1 Immunohistochemistry(IHC) 15420