

PIWIL1 Polyclonal Antibody



Catalog Number: 161042

Product name

PIWIL1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to PIWIL1

Preparation

Antigen: Recombinant protein of mouse PIWIL1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - PIWIL1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PIWIL1 antibody . Blue: DAPI for nuclear staining.