

RAB8A Polyclonal Antibody



Catalog Number: 163508

Product name

RAB8A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RAB8A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human RAB8A (NP_005361.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - RAB8A Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using RAB8A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.