

Anti-SNF2H Rabbit antibody



Catalog Number: 166826

Product name

Anti-SNF2H Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SNF2H Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SNF2H

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000

ICC/IF: 1/50

Validations



Western blot detection of SNF2H in K562, C6, 3T3, HeLa cell lysates using SNF2H Rabbit pAb(1:1000 diluted). Predicted band size: 122kDa. Observed band size: 122kDa.