

Anti-Peroxiredoxin 5 Rabbit antibody



Catalog Number: 166852

Product name

Rabbit IgG

Anti-Peroxiredoxin 5 Rabbit antibody

Applications

Specificity

WB, IHC-P, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

Peroxiredoxin 5 Rabbit polyclonal antibody

IHC: 1/50

Preparation

FC: 1/20

Antigen: A synthetic peptide of human Peroxiredoxin 5

IP: 1/20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type



Western blot detection of Peroxiredoxin 5 in HeLa cell lysates using Peroxiredoxin 5 Rabbit pAb(1:1000 diluted). Predicted band size:22kDa.Observed band size:17kDa.