

Anti-Collagen I Rabbit antibody



Catalog Number: 166894

Product name

Anti-Collagen I Rabbit antibody

Specificity

Human, Cow

Antibody description

Collagen I Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Collagen I

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

Validations



Western blot analysis of Collagen I in Human stomach tissue lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using Collagen I Antibody.