

Anti-Phospho-NF-kB p65 (Ser529) Rabbit antibody



Catalog Number: 167123

Product name

Anti-Phospho-NF-kB p65 (Ser529) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-NF-kB p65 (Ser529) Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around phosphorylation site of serine 529 (L-L-S(p)-G-D) derived from Human NFκB-p65.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

IgG


Applications

WB, IHC-P, ICC/IF


Dilutions

WB:1/500-1/1000 ICC:1/100-1/200 IHC:1/50-1/100


Validations



Immunofluorescence staining of methanol-fixed HeLa cells using NFκB-p65 (Phospho-Ser529) Antibody .



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NFκB-p65 (Phospho-Ser529) Antibody (left) or the same antibody preincubated with blocking peptide (right).




peptide (right).

Western blot detection of NFκB-p65(Phospho-Ser529) in 3T3 and HeLa cells untreated or treated with Alkaline phosphatase using NFκB-p65(Phospho-Ser529) Rabbit pAb (dilution 1:500, upper) or u03b2-Actin Mouse mAb (200068-8F10, lower). Predicted band size: 65kDa. Observed band

Anti-Phospho-NF-kB p65 (Ser529) Rabbit antibody




Catalog Number: 167123



size:65kDa.

Western blot analysis of extracts from MCF cells, treated with TNF α +CA or calf intestinal phosphatase (CIP), using NF κ B-p65 (Phospho-Ser529) Antibody .

(Phospho-Ser529) Antibody .



Western blot analysis of extracts from 293 cells untreated or treated with TNF- α using NF κ B-p65 (Phospho-Ser529) Antibody .