

# Anti-Phospho-FAK (Tyr861) Rabbit antibody



Catalog Number: 167135

## Product name

WB

Anti-Phospho-FAK (Tyr861) Rabbit antibody

## Dilutions

## Specificity

WB:1/500-1/1000

Human, Mouse, Rat

## Antibody description

Phospho-FAK (Tyr861) Rabbit polyclonal antibody

## Preparation

Antigen: Peptide sequence around phosphorylation site of tyrosine 861 (H-I-Y(p)-Q-P) derived from Human FAK.

## Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Polyclonal

## Ig Type


IgG

## Applications

## Validations



Western blot analysis of extracts from A549 cells, untreated or treated with EGF(10ng/ml,10min), using FAK (Phospho-Tyr861) Rabbit pAb (167135,1:500 diluted,upper) and FAK Mouse mAb (200899,lower).



Western blot analysis of extracts from HT29 cells untreated or treated with EGF using FAK (Phospho-Tyr861) Antibody .