

# Anti-P70 S6 Kinase beta/SRK Rabbit antibody



Catalog Number: 167396

## Product name

Anti-P70 S6 Kinase beta/SRK Rabbit antibody

WB: 1/2000-1/10000

IHC: 1/20

## Specificity

Human, Mouse, Rat

ICC/IF: 1/50

FC: 1/20

## Antibody description

P70 S6 Kinase beta/SRK Rabbit polyclonal antibody

## Validations

## Preparation

Antigen: Recombinant protein of human P70 S6 Kinase beta

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions



Immunohistochemistry of P70 S6 Kinase beta in paraffin-embedded Human colon cancer tissue using P70 S6 Kinase beta Rabbit pAb at dilution



Western blot detection of P70 S6 Kinase beta in Rat Brain,C6,HeLa cell lysates using P70 S6 Kinase beta Rabbit pAb(1:1000 diluted).Predicted band size:54kDa.Observed band size:54kDa.