

Anti-LC3B Rabbit antibody



Catalog Number: 167419

Product name

Anti-LC3B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LC3B Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human LC3B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/100

IP: 1/50

Validations



Western blot detection of LC3B in U87-MG, HeLa, C6 and 3T3 cell lysates using LC3B Rabbit pAb (1:1000 diluted). Predicted band size: 15kDa. Observed band size:14, 16kDa.