

Anti-GRB2 Rabbit antibody



Catalog Number: 167501

Product name

Anti-GRB2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

GRB2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human GRB2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/2000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of GRB2 in MCF7,3T3,Raw264.7,HeLa cell lysates using GRB2 Rabbit pAb(1:1000 diluted).Predicted band size:25KDa.Observed band size:25KDa.