

Anti-GRB2 Rabbit antibody



Catalog Number: 167501

Product name

Rabbit IgG

Anti-GRB2 Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000

GRB2 Rabbit polyclonal antibody

IHC: 1/50

Preparation

IP: 1/20

Antigen: A synthetic peptide of human GRB2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of GRB2 in MCF7,3T3,Raw264.7,Hela cell lysates using GRB2 Rabbit pAb(1:1000 diluted).Predicted band size:25KDa.Observed band size:25KDa.