

Anti-P4HB Rabbit antibody



Catalog Number: 167502

Product name

Anti-P4HB Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

P4HB Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human P4HB

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/20

IP: 1/20

Validations



Western blot detection of P4HB in K562,C6 cell lysates using P4HB Rabbit pAb(1:1000 diluted).Predicted band size:57kDa.Observed band size:57kDa.