

# Anti-BRK1 Rabbit antibody



Catalog Number: 167526

## Product name

Anti-BRK1 Rabbit antibody

## Specificity

Human

## Antibody description

BRK1 Rabbit polyclonal antibody

## Preparation

Antigen: Synthetic peptide of human BRK1

## Formulation

-20°C, pH7.4 PBS, 0.05% NaN<sub>3</sub>, 40% Glycerol

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC-P

## Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

## Validations



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 167526(BRK1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 167526(BRK1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original



magnification: ×200)

Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane: A172 cells, Primary antibody: 167526(BRK1 Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1

# Anti-BRK1 Rabbit antibody



Catalog Number: 167526

---

minute