

Anti-Histone H3.3 Rabbit antibody



Catalog Number: 167531

Product name

Anti-Histone H3.3 Rabbit antibody

Dilutions

WB: 1/1000

Specificity


Human, Mouse, Rat

IHC: 1/200

Antibody description

Histone H3.3 Rabbit polyclonal antibody

Validations



Preparation

Antigen: A synthetic peptide of human Histone H3.3

Western blot detection of Histone H3.3 in K562, C6, 3T3 cell lysates using Histone H3.3 Rabbit pAb(1:1000 diluted). Predicted band size: 15kDa. Observed band size: 15kDa.

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P



Immunohistochemistry of Histone H3.3 in paraffin-embedded Human lung cancer tissue using Histone H3.3 Rabbit pAb at dilution 1/50