

Anti-Calcineurin A Rabbit antibody



Catalog Number: 167566

Product name

Anti-Calcineurin A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Calcineurin A Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Calcineurin A

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, FC, IP


Dilutions

WB: 1/2000

FC: 1/20

IP: 1/20

Validations



Western blot detection of Calcineurin A in HeLa, A549, U2OS cell lysates using Calcineurin A Rabbit mAb(1:1000 diluted). Predicted band size: 59kDa. Observed band size: 59kDa.