

# Anti-Daxx Rabbit antibody



Catalog Number: 167649

## Product name

Anti-Daxx Rabbit antibody

## Specificity

Human

## Antibody description

Daxx Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Daxx

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/2000

IHC: 1/20

ICC/IF: 1/20

FC: 1/20

## Validations



Western blot detection of Daxx in  
HepG2,U251,HL-60 using Daxx Rabbit  
mAb(1:1000 diluted)